

The MASTERS' Collection gives you the freedom to select unique finishes, materials, colors, and designs – while staying within your budget. The most expansive array of materials and finishes ever seen in the market place. Virtually unlimited combinations to create your masterpiece.

ICON MASTERS' COLLECTION

The ICON makes a bold aesthetic statement in any treatment room. Add our optional comfort features to create an aesthetic and sensory experience that will not be forgotten.



SPECIFICATIONS

- Removable and replaceable padding
- AeroCel™ padding for ultimate comfort
- Large height range - 25" - 35"



TABLE SPECIFICATIONS

Lift System	Electric
Weight Capacity	500 lbs. (226.5 kg)
Fabric	TerraTouch or Ultraleather
Finish	Laminate
Padding	4" (10.2 cm) or 6" (15.2 cm) AeroCel™ (Optional)
Height Range 4" (10.2 cm) pad	25" - 35" (62 cm - 88 cm)
Height Range 6" (15.2 cm) pad	27" - 37" (67 cm - 93 cm)
Width	31" (79 cm), 33" (84 cm)
Length	75" (190.5 cm)
Options (must be ordered with table)	SoVi Vibroacoustic System, Fuzion Heat, ABC System, LED Lighting, External Outlet
Warranty	3 years parts, 2 years labor
Safety Listings	cETLus listed 3034177

Accessories

Adjustable Side Arm Rests, Boiance Shelf Support, Boiance Face Rest Cushion, Quick Lock Face Rest

LAMINATE COLORS



FABRIC COLORS

Meets the requirements of ISO 10993-1:2009 Biological Evaluation of Medical Tables & passes testing for Cytotoxicity, Skin Irritation & Sensitization as well as CA TB117 Fire Retardance testing.

To view our current list of available fabric colors, scan QR code or visit spatables.com



OPTIONS

ABC SYSTEM


Oakworks Patented ABC system allows you to create the perfect level of support for each client, taking pressure off their breasts and tension out of their back. The system also raises the pad above level, making it a great head support for face up positions, and can be used for stretches also.



Incredible breast comfort
Lower to lessen pressure




Electric adjustments
Large adjustment range



SO•VI

Optional SoVi vibroacoustic system is designed to deliver all the latest forms of healing modalities via a very flexible integrated delivery system.



FUZION HEAT

Optional Fuzion heating system is embedded in the foam, creating a more even heat.



717.235-6807 www.spatables.com email: customerservice@oakworks.com