



# Procedure Chairs

3000 | 3050 | 3100



## Mobility System

Easily move your chair in the room  
(not for patient transport)



# 3050 Procedure Chair



**3050 Model**  
Arm Rest and Mobility  
System standard  
(not for patient transport)


- Powered: height, back rest  
leg rest, tilt, zero gravity
- ADA 19" low height
- 550 lb weight rating



Oakworks Procedure Chairs, with their 19" low height and fully powered functions are a cost effective ergonomic solution for many procedures and exams. Built in our Pennsylvania factory and backed by a 3 year parts, 2 years labor warranty.



**3000 Model**  
(shown with optional  
headrest)



Arm rests can be rotated horizontally and can be folded down, out of the way, providing excellent access



Our Mobility System at the push of a button allows you to move the chair and lock it in place.



**see page 4 for specifications on each model**

# 3100 Procedure Chair

- Powered: height, back rest, leg rest, tilt, zero gravity
- ADA 19" low height
- 550 lb weight rating



## 3100 Model

Mobility System standard  
(not for patient transport)

The 3100 Procedure Chair comes standard with stirrups and a narrower leg section that permits the stirrups to be extended or retracted with the patients legs supported.

### Flat Headrest Pad



Adjustable in height and angle, U-Shaped Pad for prone available

### Hand Control




User programmable memory positions and presets

### Removable Pads



Remove for cleaning or replacement

### Seat Section



Our unique hinge design shortens the seat for more comfort.


see page 4 for specifications on each model

# Comparison Chart/ Specifications



All 3000 Series Chairs can be configured to your exact needs with many options and accessories available. They come fully assembled with a 3 year parts, 2 year labor warranty.

O\* - Arm Rest or Side Rails - not both



Spec	Description	3000	3050	3100
3 Section Top	26" W x 61" L (76" w/head rest)	S	S	S
Backrest	25" L	S	S	S
Seat	24" L	S	S	S
Leg Rest	11" L x 26" W (19" W on 3100)	S	S	S
Motion 1	Height Range 19" - 33"	S	S	S
Motion 2	Trendelenburg 15° / Reverse 21°	S	S	S
Motion 3	Backrest 0-75°	S	S	S
Motion 4	Leg Rest 0-90°	S	S	S
Weight Capacity	550 lb	S	S	S
Hand Control	Programmable Memory	S	S	S
Foot Control	All functions	O	O	S
Padding	2.5" Comfort Foam™	S	S	S
Flat Headrest	11" L x 13" W (extends up to 4")	O	S	S
Stirrups	Adjustable / Retractable	N/A	N/A	S
Procedure Tray	12.25" x 10.25" (holds 2 quarts)	N/A	N/A	O
Mobility System	Electronic Caster System	N/A	S	S
Side Rails	25.5" L x 6" Tall (pair)	O*	N/A	O*
Arm Rests	4" wide x 18" L (pair)	O*	S	O*
Battery Backup	Automatic backup	O	O	O
Paper Roll Holder	Accommodates rolls up to 18" W	S	S	S
Safety Listings	FDA listed, CE marked, cETLus listed, CB Certificate	S	S	S



**View  
Product  
Video**

**small  
business**  
Registration #130256837

**FDA LISTED  
CLASS I PRODUCT**



## FABRIC COLORS

White-T26	Coal-T01
Opal-T12	Espresso-T27
Heron-T05	Ruby-T08
Earth-T07	Ocean-T20
Stone-T13	Pewter-T29
Sapphire-T17	Sage-T22
Clay-T18	

www.oakworksmed.com 717-235-6807

© 2023 OAKWORKS®, Inc. 923 East Wellspring Rd. New Freedom, PA 17349 USA