




The MASTERS' Collection gives you the freedom to select unique finishes, materials, colors, and designs – while staying within your budget. The most expansive array of materials and finishes ever seen in the market place. Virtually unlimited combinations to create your masterpiece.

TALISE MASTERS' COLLECTION

The TALISE table has simple lines and curved panels for a sense of softened strength. Delivers the ultimate comfort level and ergonomics for clinicians and clients. Two height range options and three cabinet options.



SPECIFICATIONS


- ADA Compliant Low Height or High Height options
- Choose finishes for the cabinet and corners
- 3 cabinet styles - standard, storage and warming drawer

TABLE SPECIFICATIONS

| | |
|--------------------------------------|---|
| Lift System | Electric |
| Weight Capacity | 500 lbs. (226.5 kg) |
| Fabric | TerraTouch or Ultraleather |
| Finish | Laminate |
| Padding | 4" (10.2 cm) or 6" (15.2 cm) AeroCel™ (Optional) |
| Height Range 4" (10.2 cm) pad | 18.9" - 29.7" low height (48 cm - 75 cm) or 20.5" - 34" (52 cm - 86 cm) |
| Width | 31" (79 cm), 33" (84 cm) |
| Length | 75" (190.5 cm) |
| Options (must be ordered with table) | SoVi Vibroacoustic System, Fuzion Heat, ABC System, LED Lighting, External Outlet |
| Warranty | 3 years parts, 2 years labor |
| Safety Listings | cETLus listed 3034177 |

| | |
|-------------|---|
| Accessories | Adjustable Side Arm Rests, Boiance Shelf Support, Boiance Face Rest Cushion, Quick Lock Face Rest |
|-------------|---|

LAMINATE COLORS



FABRIC COLORS

Meets the requirements of ISO 10993-1:2009 Biological Evaluation of Medical Tables & passes testing for Cytotoxicity, Skin Irritation & Sensitization as well as CA TB117 Fire Retardance testing.

To view our current list of available fabric colors, scan QR code or visit spatables.com




FEATURES/OPTIONS

ABC SYSTEM


Oakworks Patented ABC system allows you to create the perfect level of support for each client, taking pressure off their breasts and tension out of their back. The system also raises the pad above level, making it a great head support for face up positions, and can be used for stretches also.



Incredible breast comfort
Lower to lessen pressure




Electric adjustments
Large adjustment range




SO•VI

Optional SoVi vibroacoustic system is designed to deliver all the latest forms of healing modalities via a very flexible integrated delivery system.



FUZION HEAT

Optional Fuzion heating system is embedded in the foam, creating a more even heat.



717.235-6807 www.spatables.com email: customerservice@oakworks.com