

One Portable Massage Package

A lightweight wood table specifically constructed for the massage professional whose practice requires constant travel.

Package Includes

- One table in natural wood finish
- Shiatsu Cables and Access End
- Adjustable Face Rest Platform
- AeroCel Face Rest Crescent
- Arm Hammock



Shown with Coal top



PATENTED
INTEGRAHINGE™



UNILOCK™



ULTIMATE ACCESS END PANEL
& SHIATSU CABLES



CABLE LOCK
SYSTEM

TERRATOUGH COLORS

Meets the requirements of ISO 10993-1:2009 Biological Evaluation of Medical Tables & passes testing for Cytotoxicity, Skin Irritation & Sensitization as well as CA TB117 Fire Retardance testing. Please refer to the product cleaning instructions for a list of approved cleaners and disinfectants.





TABLE SPECIFICATIONS

Max. Operating Weight	550 lbs.
Weight	25-32 lbs.
Height Ranges	24"-34"
Widths	29", 31"
Length	73"
Tops	Rectangular
Padding	Plush
Standard	Access End Panel, Shiatsu Cables
Accessories	QuickLock Face Rest Platform, Arm Hammock, and our full line of accessories
Fabric	TerraTouch

Contact Oakworks
717.235.6807

www.massagetables.com
customerservice@oakworks.com