

# Portal Pro Massage Chair

The patented PortalPro<sup>®</sup> has been the undisputed favorite massage chair for over 40 years. This lightweight chair offers superior ergonomics, flexibility and strength. With a working weight of 350 lbs, the PortalPro<sup>®</sup> easily adjusts to fit any size client.



STERNUM PAD OPTION



ERGONOMIC DESIGN



ACCESSIBILITY AND  
ADJUSTABILITY



WHEELS



Includes Carry Case w/Strap,  
Shown with Ocean Fabric

## TERRATOUGH COLORS

Meets the requirements of ISO 10993-1:2009 Biological Evaluation of Medical Tables & passes testing for Cytotoxicity, Skin Irritation & Sensitization as well as CA TB117 Fire Retardance testing. Please refer to the product cleaning instructions for a list of approved cleaners and disinfectants.



## TABLE SPECIFICATIONS

<i>Max. Operating Weight</i>	<i>350 lbs.</i>
<i>Weight</i>	<i>19 lbs.</i>
<i>Frame</i>	<i>Black Powder-Coated Aluminum Frame</i>
<i>Standard</i>	<i>Carry Case with Strap</i>
<i>Fabric</i>	<i>TerraTouch</i>
<i>Accessories</i>	<i>Sternum Pad</i>

Contact Oakworks  
717.235.6807  
[www.massagetables.com](http://www.massagetables.com)  
[customerservice@oakworks.com](mailto:customerservice@oakworks.com)