



POWERLINE EXAM TABLES



Sapphire fabric, Grey finish
and optional storage shelf



Sapphire fabric, Natural Finish
and optional storage shelf

Visit www.oakworksmed.com for a quote!

Features

- Cost effective solution for any facility
- Partially assembled for quick set up
- Hardwood construction has a 500 lb. weight capacity
- 3 year parts warranty
- Available in 24", 27" and 30" widths
- Flat top and backrest options
- Highly water resistant and can be wiped down with disinfectants as needed

WAVE BACKREST



The wavedesign makes adjusting your backrest quick and easy.

DURABLE MATERIALS



Hardwood construction, reinforced hardware, durable padding & fabric.


PAPER ROLL HOLDER



Paper Roll Holder is included standard. Paper roll not included.

TERRATOUGH COLORS

Meets the requirements of ISO 10993-1:2009 Biological Evaluation of Medical Tables & passes testing.



FINISHES



WEDGE BOLSTER

This product is not FDA listed.

Easily used as a backrest or leg support. 24"W x 10"H x 24"L with a 36° incline.



CONFIGURE YOUR TABLE

Backrest or Flat Top	24" (61cm), 27" (69cm) or 30" (76cm) wide x 72" (183cm) long
Height	30" (76cm)
Backrest Top	Manual Multi-Position
Weight Capacity	500 lb (227kg)
Product Weight	87 lb (39kg)
Padding	1.75" (4cm) Firm Response or 2.5" (6cm) Comfort Foam™
**Base Finish	Choice of: Natural and Grey
Paper Roll Holder	fits up to 22" (56cm) wide paper rolls
Storage Shelf	Fits all widths
Warranty	3 years parts, labor not included
Safety Listings	FDA listed

ACCESSORIES order any time.		#
PowerLineTableExtender	27" (69cm) wide x 6" (15cm) long	51909-T*
PowerLineTableExtender	30" (76cm) wide x 6" (15cm) long	51910-T*
QuickLock™ Head Rest & AeroCel Pad	Comfortable, adjustable prone or supine support	PKG7554-T*
Stabilization Strap	Stabilize your patient	63530
Prone Pillow	For prone or supine positioning	66180-T01

*Fabric color code must be specified.

** Finish color must be specified.

PRONE PILLOW

- Increased ventilation & improved comfort
- Low profile reduces cervical extension/flexion
- Superior positioning for any size patient
- Removable, easy to clean pads



717.759.3146 www.oakworksmed.com email: medsales@oakworks.com