


Wellspring

Our lightest, full size table features our patented complementary support system and the best access on the market. This stylish lightweight table is supported by a sturdy ultimate access aluminum understructure ensuring superior leg room for the therapist on all sides of the table. This is a perfect table for anyone who travels, treats a diverse client base, or works seated much of the time.



Shown in Sapphire



SUSPENSION SYSTEM



UNILOCK™



AERO • CEL™ PADDING UPGRADE



TELESCOPING LEGS




EXCEPTIONAL FRAME STRENGTH




PATENTED INTEGRATING™

TERRATOUGH COLORS

Meets the requirements of ISO 10993-1:2009 Biological Evaluation of Medical Tables & passes testing for Cytotoxicity, Skin Irritation & Sensitization as well as CA TB117 Fire Retardance testing. Please refer to the product cleaning instructions for a list of approved cleaners and disinfectants.



Heron



Ocean



Sapphire




Espresso



Pewter



Coal




Opal



Clay



Sky Blue



Stone



White

TABLE SPECIFICATIONS

Max. Operating Weight	550 lbs.
Weight	25-30 lbs.
Height Ranges	24"-34"
Widths	29", 31"
Length	73"
Tops	Rectangular, Round Corner
Padding	Plush, AeroCel, Semi-Firm
Standard	Ultimate Access Panel
Accessories	QuickLock Face Rest Platform, Arm Hammock, and our full line of accessories
Fabric	TerraTouch

Contact Oakworks
717.235.6807
www.massagetables.com
customerservice@oakworks.com