

OAKWORKS®

THE POWER OF OAKWORKS®



THE POWER OF STRENGTH

made in the USA
with US & imported parts



ATHLETIC TRAINING EQUIPMENT

717.759.3146

www.oakworksmed.com

TABLE OF CONTENTS



PAGE 1

PowerLine™ Treatment Table

PAGE 2

PT100 Treatment Table

PAGE 3

BOSS™ Portable Treatment Table

PAGE 4

Portable Manipulation Table

PAGE 5

Storable Mat

PAGE 6

Portable Taping Table

PAGE 7

Prone Pillow

PAGE 7

Bolsters

PAGE 8

Stools

PAGE 8

Face Rests and covers

PAGE 9

Fabric

BACK COVER

Featured Product - Stabilization Straps

made in the USA
with US & imported parts



SAFETY AND TESTING

We know that your clients, your career and your livelihood depend on the safety, reliability and quality of the products you use. Simply put, nobody else tests like OAKWORKS®, Inc. Every one of our new designs goes through a rigorous, specialized, four-part testing protocol including a dynamic load test, a lateral force test, a stability test and component and cycle tests. This equipment is engineered to exceed all the challenges your table might experience in the field, before it gets there.



Dynamic load testing simulates real-life hop and drop table use.

OAKWORKS®

POWERLINE™



Now available
in three widths!
24", 27" and 30"



500 LB.
BARIATRIC
WEIGHT
CAPACITY

Three widths available - 24", 27", 30"

Padding - Comfortable 2.5" padding

Colors - 11 fabric colors available

Included feature - Paper roll holder

Rock solid all hardwood base construction with no particle board; built to last.

Ships partially assembled & takes minutes to set up - Tools Included!


Strongest table in its class with a Bariatric Rated 500 lb Weight Capacity!



Finishes



GREY



NATURAL

Highly water resistant finish designed to hold up to years of disinfecting and cleaning.

made in the USA
with US & imported parts



3 YEAR WARRANTY!

Specifications

| | |
|-------------|--------------------------------------------------------------------------------------------|
| LOAD WEIGHT | Bariatric Rated 500 lbs |
| HEIGHT | 30" |
| WIDTH | 24", 27" or 30" |
| LENGTH | 72" |
| TOP | Flat Top or Backrest |
| FABRIC | 11 TerraTouch™ Colors |
| PADDING | 2.5" Semi-Firm |
| STANDARD | Accessory Outlets, Paper Roll Holder |
| ACCESSORIES | QuickLock™ Face Rest Platform, Aero-Cel™ Face Rest Crescent, Table Extender, Storage Shelf |

Italic = Additional Charge

PT100

Super Stiff Design is more stable than the competition!

PT100



16"-34"
Height Range!



Heavy Gauge Welded Steel Dual Crossbeam construction lifts the table vertically, taking up less room than cantilevered designs.

Optional Locking Casters ensure safety and stability.



Accessory Outlets allow versatility for accessories.

made in the USA
with US & imported parts



Specifications

| | |
|-------------|-------------------------------------------------------------------------|
| LOAD WEIGHT | 550 lbs |
| HEIGHT | 16" - 34" (with 1.75" foam option) |
| WIDTH | 27" or 29" |
| LENGTH | 73" |
| TOP | Flat Top |
| FABRIC | 11 TerraTouch™ Colors |
| PADDING | 1.75" Firm Response, 2.5" Comfort Foam |
| CONTROLS | Foot Control |
| REGULATORY | FDA listed, ETL listed |
| STANDARD | Accessory Outlets |
| OPTIONS | <i>Locking Twin Casters (one end), Hand Control, Extra Foot Control</i> |

Italic = Additional Charge

BOSS™ *Take advantage of the Golden Minute!*

BOSS



made in the USA
with US & imported parts



All weather design holds up
on the sideline for many years.

Field Feet - Removable feet
allow a firm footing on any surface.



*The **STRONGEST** and most **DURABLE**
portable table on the market –*

THE BOSS™ PACKAGE

Includes carry case and field feet



Specifications

| | |
|--------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| LOAD WEIGHT | 500 lbs |
| TABLE WEIGHT | 32-35 lbs |
| HEIGHT | 24" - 34" |
| WIDTH | 29" or 31" |
| LENGTH | 73" |
| TOP | Flat Top |
| FABRIC | Tufflex™ Black or Royal Blue |
| PADDING | 2.5" Semi-Firm |
| STANDARD | Patented Intergra-Hinge™, UniLock™, Accessory Outlets, Complementary Suspension System, Exceptional Frame Strength, Telescoping Aluminum Legs, Field Feet |
| OPTIONS | 11 TerraTouch™ Colors |
| ACCESSORIES | QuickLock™ Face Rest Platform, Aero-Cel™ Face Rest Crescent, Carry Case |

Italic = Additional Charge

CS System™ (Complementary Suspension) : U.S. Patent #6,192,809

UniLock™ : U.S. Patent # 7,472,656

See the BOSS™ in action...

ESPN X-Games • Memphis Grizzlies • Nashville Predators • Atlanta Braves • Charlotte Bobcats • Detroit Tigers • Disney Animal Kingdom Show Support • Barry University • Bowie State University • Campbellsville University • Emory University • Grand Canyon University • Harvard University • Indiana State University • Penn State • Rowan University • Temple University • Kennesaw State University • Post University • Rider University-Alumni Gym • St. Joseph's College • Stevenson University • University of Findlay • University of Mass (Dartmouth) • University of Oregon • University of Wyoming • Villanova University Athletics • West Virginia University of Tech • Wofford College ... AND MANY MORE!

PORTABLE MANIPULATION TABLE - METAL



made in the USA
with US & imported parts



Light weight and compact - the Portable Manipulation table is the choice when on the go. The narrower design supports side stretches and great access.

Specifications

| | |
|-----------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|
| LOAD WEIGHT | 500 lbs |
| TABLE WEIGHT | 26 lbs |
| HEIGHT | 24" - 34" |
| WIDTH | 27" |
| LENGTH | 73" |
| TOP | Flat Top |
| FABRIC | 11 TerraTouch™ Colors |
| PADDING | 1.5" Firm Response |
| STANDARD | Patented Intergra-Hinge™, UniLock™, Accessory Outlets, Complementary Suspension System, Exceptional Frame Strength, Telescoping Aluminum Legs |
| ACCESSORIES | <i>QuickLock™ Face Rest Platform, Aero-Cel™ Face Rest Crescent, Carry Case</i> |
| <i>Italic = Additional Charge</i> | |

UniLock™: U.S. Patent # 7,472,656



"The Portable Manipulation Table from Oakworks offers an outstanding platform for performing thrust & non-thrust manipulation. The team at Oakworks has listened to the feedback of practicing manual physical therapists & created a high value & excellent performing product. You can't ask anything more out of a portable manipulation table!"

Timothy W. Flynn, PT, PhD



ADJUSTABLE

The 10" height range supports good ergonomics for everyone.



CARRY CASE

Take your table with you to home visits & maintain a consistent working surface for measurable results.

STORABLE MAT

STORABLE MAT



Ultra Stiff Leg Bracing holds up under heavy lateral force.

Space saving design folds up easily for storage and quick portability.

Stabilizing Feet for maximum stability. (standard on 40" width)

Choose your width:

33", 35" or 40"

The ONLY portable mat table in the market with a 550 lb Weight Capacity!



made in the USA
with US & imported parts



Specifications

| | |
|-----------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|
| LOAD WEIGHT | 550 lbs |
| TABLE WEIGHT | 40-55 lbs (depending on width) |
| HEIGHT | 16" - 25" |
| WIDTH | 33", 35", 40" |
| LENGTH | 78" |
| TOP | Flat Top |
| FABRIC | 11 TerraTouch™ Colors |
| PADDING | 2.5" Semi-Firm, 1.5" Firm |
| STANDARD | Patented Intergra-Hinge™, UniLock™, Accessory Outlets, Complementary Suspension System, Exceptional Frame Strength, Telescoping Aluminum Legs |
| OPTIONS | Stabilizing Feet for 33" & 35" (standard on 40" width) |
| ACCESSORIES | QuickLock™ Face Rest Platform, Aero-Cel™ Face Rest Crescent, Carry Case (not available for 40" width) |
| <i>Italic = Additional Charge</i> | |

CS System™ (Complementary Suspension) : U.S. Patent #6,192,809
UniLock™: U.S. Patent # 7,472,656

PORTABLE TAPING TABLE

PORTABLE TAPING TABLE

Marine-Grade Polymer Top - super tough and handles heavy loads with ease.

Strong, sturdy design offers great sideline utility and durability.

Portability - Easy to store and carry with convenient folding design.

Field Feet - Removable feet allow a firm footing on any surface.



*Shoulder Strap
Included!*



made in the USA
with US & imported parts



*The world's BEST portable adjustable
height taping table! 32" - 42" Height Range!*

THE PTT PACKAGE

Includes carry case and field feet

Specifications

| | |
|--------------|-------------------------------------------------------|
| LOAD WEIGHT | 500 lbs |
| TABLE WEIGHT | 35 lbs. |
| HEIGHT | 32" - 42" |
| WIDTH | 24" |
| LENGTH | 48" |
| TOP | Marine-Grade Polymer Top |
| STANDARD | Telescoping Aluminum Legs, Field Feet, Shoulder Strap |
| ACCESSORIES | <i>Carry Case</i> |

Italic = Additional Charge



PRONE PILLOW & HYGIENIC COVERS



available in TerraTouch™ Coal only



- Superior positioning for any size/shape
- Large vent areas for easy breathing
- Removable pads are easy to clean & replaceable
- Most comfortable design on the market

PRONE PILLOW

BOLSTERS

100% Polyurethane fabric, no PVC, handle on the end

Available in 11 TerraTouch™ Colors!



★ **most popular** ★



WEDGE BOLSTER

24"W x 10"H x 24"L with a 36° Incline.

made in the USA
with US & imported parts



BOLSTERS

For our full line of bolsters, please visit www.oakworksmmed.com

STOOLS

OAKWORKS® Stools are durable, comfortable and support 300 lbs. with simple pneumatic height adjustments.

BASIC, PREMIUM & PROFESSIONAL STOOLS

STOOLS



BASIC STOOL

Colors: Coal ONLY
Height Range: 18" - 24"



PREMIUM STOOL

Colors: All TerraTouch™ Colors
Low Height: 18" - 24"
High Height: 20" - 28"



PREMIUM STOOL with SADDLE SEAT

Colors: All TerraTouch™ Colors
Low Height: 19" - 25"
High Height: 22" - 30"



PREMIUM STOOL with BACKREST

Colors: All TerraTouch™ Colors
Low Height: 18" - 24"
High Height: 20" - 28"

ACCESSORIES

QUICKLOCK™ FACE REST PLATFORM

Rock solid support over a wide range of positions via double cam lock system.



AERO-CEL™ FACE REST CRESCENT

Soft and supportive, this versatile face rest crescent adjusts to support any size or shape face for prone comfort.



FLAT DISPOSABLE FACE REST COVERS

100 per box. Super soft and disposable, this flat cover provides a sanitary covering.



FITTED DISPOSABLE FACE REST COVERS


50 per box. Provides a fitted disposable sanitary covering at an irresistible price.



FABRIC

TerraTouch™

Meets the requirements of ISO 10993-1:2009 Biological Evaluation of Medical Tables & passes testing for Cytotoxicity, Skin Irritation & Sensitization as well as CA TB117 Fire Retardance testing.



White



Stone



Clay



Ocean



Pewter




Opal



Heron




Sapphire



Coal



Sky Blue



Espresso

Tufflex™

Specifically designed for use on the sidelines, Tufflex can withstand cleats and sun.

Black or Royal Blue are standard on the BOSS™ table.



ROYAL BLUE



BLACK

OAKWORKS®

STABILIZATION STRAPS




- Adjustable in length over a wide range for stabilization of any body part
- Designed to work with OAKWORKS® PowerLine™ stationary tables and all of OAKWORKS® portable tables except the Storable Mat and Portable Taping Table
- Secure and strong, heavy duty design will last for years
- Wide straps offer more surface area for client comfort
- Great for range of motion, mobilization, isolation



Watch the Video!

featuring Dr. Robert DuVall
PT, DHSc, ATC, OCS, FAAOMPT CSCS



**THE POWER
OF STRENGTH**

OAKWORKS®