



Oakworks Imaging Tables • Motions At-A-Glance



Trendelenburg Tilt



Longitudinal Travel



Lateral Tilt



Lateral Travel

Table Type	Height Range *without pad	Weight Rating	Trendelenburg Tilt	Longitudinal Travel	Lateral Tilt	Lateral Travel
CFPM100	22"-40" (56-102cm)	500lbs	n/a	n/a	n/a	n/a
CFPM200	25"-43" (64-109cm)	500lbs	n/a	n/a	±15°	n/a
CFPM201	25"-42" (64-107cm)	500lbs	+15°/-11°	n/a	n/a	n/a
CFPM300	26"-43" (66-109cm)	500lbs	+15°/-11°	n/a	±15°	n/a
CFPM301	25"-43" (64-109cm)	500lbs	+15°/-11°	10" (25cm)	n/a	n/a
CFPMB301	26-44" (66-112cm)	750lbs	+15°/-11°	10" (25cm)	n/a	n/a
CFPM302	25"-43" (64-109cm)	500lbs	n/a	10" (25cm)	±15°	n/a
CFPM400	26"-43" (66-109cm)	500lbs	+15°/-11°	10" (25cm)	±15°	n/a
CFPM401	26"-44" (66-112cm)	500lbs	+15°/-11°	10" (25cm)	n/a	8" (20cm) ±4"
DTPM300	24"-40" (61-102cm)	450lbs	+14°/-14°	23" (58cm) (manual)	n/a	n/a

Table Top Specifications - 56" Metal Free Imaging Space

