



Procedure Chairs

3000 | 3050 | 3100

made in the USA
with US & imported parts



Mobility System

Easily move your chair in the room
(not for patient transport)

19"



*A cost-effective, ADA
compliant procedure chair.*

OAKWORKS®
MEDICAL PRODUCTS

3050 Procedure Chair

made in the USA
with US & imported parts



3050 Model
Arm Rest and Mobility
System standard
(not for patient transport)

- Powered: height, back rest, leg rest, trendelenburg
- ADA 19" low height*
- 550 lb weight rating


* high height: 22" - 36" (56 - 91 cm) also available



3000 Model
(shown with optional
headrest)

Oakworks Procedure Chairs, with their 19" low height* and fully powered functions are a cost effective ergonomic solution for many procedures and exams. Built in our Pennsylvania factory and backed by a 3 year parts, 2 years labor warranty.

* high height: 22" - 36" (56 - 91 cm) also available



Arm rests can be rotated horizontally and can be folded down, out of the way, providing excellent access



Our Mobility System at the push of a button allows you to move the chair and lock it in place.



see page 4 for specifications on each model

3100 Procedure Chair

- Powered: height, back rest, leg rest, trendelenburg
- ADA 19" low height*
- 550 lb weight rating

* high height: 22" - 36" (56-91 cm) also available


made in the USA
with US & imported parts



3100 Model

Mobility System standard
(not for patient transport)


Better Access, Better Ergonomics



Eliminating the large central column many use, our twin lifting tower design supports access from all sides of the chair.

The 3100 Procedure Chair comes standard with stirrups and a narrower leg section that permits the stirrups to be extended or retracted with the patients legs supported.

Flat Headrest Pad



Adjustable in height and angle, U-shaped pad for prone available

Hand Control




User programmable memory positions and presets

Removable Pads



Remove for cleaning or replacement

Seat Section



Our unique hinge design shortens the seat for more comfort

Accessory T-Rail



Easily attached to multiple locations

Side Rails



Fold down completely out of the way

Stirrups



Retractable, adjustable in length and width

Procedure Tray



Retractable, stainless steel

see page 4 for specifications on each model

Comparison Chart/ Specifications



All 3000 Series Chairs can be configured to your exact needs with many options and accessories available. They come fully assembled with a 3 year parts, 2 year labor warranty.

O* - Arm Rest or Side Rails - not both



Spec	Description	3000	3050	3100
3 Section Top	26" W x 61" L (76" w/head rest)	S	S	S
Backrest	25" L	S	S	S
Seat	24" L	S	S	S
Leg Rest	11" L x 26" W (19" W on 3100)	S	S	S
Motion 1	Height Range: 19" - 33" (48cm-84cm) low height or 22" - 36" (56-91cm) high height	S	S	S
Motion 2	Trendelenburg 15° / Reverse 21°	S	S	S
Motion 3	Backrest 0-75°	S	S	S
Motion 4	Leg Rest 0-90°	S	S	S
Weight Capacity	550 lb	S	S	S
Hand Control	Programmable Memory	S	S	S
Foot Control	All functions	O	O	S
Padding	2.5" Comfort Foam™	S	S	S
Flat Headrest	11"L x 13" W (extends up to 4")	O	S	S
Stirrups	Adjustable / Retractable	N/A	N/A	S
Procedure Tray	12.25" x 10.25" (holds 2 quarts)	N/A	N/A	O
Mobility System	Electronic Caster System	N/A	S	S
Side Rails	25.5" L x 6" Tall (pair)	O*	N/A	O*
Arm Rests	4" wide x 18" L (pair)	O*	S	O*
Battery Backup	Automatic backup	O	O	O
Paper Roll Holder	Accommodates rolls up to 18" W	S	S	S
Safety Listings	FDA listed, CE marked, cETLus listed, CB Certificate	S	S	S



**View
Product
Video**

**small
business**
Registration #130256837



**FDA LISTED
CLASS I PRODUCT**

FABRIC COLORS



www.oakworksmed.com 717-759-3146

© 2024 OAKWORKS®, Inc. 923 East Wellspring Rd. New Freedom, PA 17349 USA

Med-MKG-MKT-139 Rev.3